

## Updating Chrombox O (Optimizer)

1. Unzip the archive with the new code.
2. The folder containing the code, e.g. 'O-13-09' should be placed in the folder "...\\OO\\code".
3. Open the file "oo\_localsettings.sdv" that is found in the O root folder and update the version to the folder name of the code. The part to be edited is shown in blue in the example below.
4. An alternative to (3) is to start Chrombox O, select [Settings] down in the right corner of the window followed by [Directories]. Choose the new version under "code version" and press [Save local settings]. The new code will be active the next time the program is started.
5. Note that it is not necessary to delete the folders with old code. Keeping these allow you to run previous versions if necessary.

An example of "oo\_localsettings.sdv" is shown below. The field to check or edit is shown in blue.

---

```
defaultfolders; 1; 1 use default folder settings, 0 use paths
specified below
defaultmethod; Default; Method to load on startup
path_designs; k:\\CHROMBOX\\OO\\designs; Folder for saved designs
path_experiments; k:\\CHROMBOX\\OO\\experiments; Folder for
experimental data
path_export; k:\\CHROMBOX\\OO\\export; Folder for import/export of
various data
path_method; k:\\CHROMBOX\\OO\\methods; Folder for methods
path_models; k:\\CHROMBOX\\OO\\models; Folder for saved models
path_rawdata; k:\\CHROMBOX\\OO\\rawdata; Rawdata folder
path_reports; k:\\CHROMBOX\\OO\\reports; Folder for reports
path_results; k:\\CHROMBOX\\OO\\results; Folder for results
tracker; 0; For development purposes, 0 or 1
user; Anonymous; User ID for info fields
version; O-13-09; Code version to use
windowpos; [0.05 0.05 0.9 0.9]; Window position and size
[leftposition lowerposition width height] in fractions of
screen size
```

---